

JUST PUSH

FOR HIGH SENSITIVITY,
HIGH SPECIFICITY,
HIGH REPRODUCIBILITY,
WITH VIRTUALLY NO PAIN,
AND LOW TECHNIQUE
VARIATION BETWEEN
USERS.

UniTest® PC (pain control)

THE NEUROLOGICAL GATE-CONTROL THEORY involves activating nerve cells that block pain transmission. The UniTest® PC touch post activates these nerve cells to suppress pain, while the refined points provide well-defined wheals with very little tissue trauma.

Allow touch post to contact skin for one second then press firmly into skin. ▶

The UniTest® PC TOUCH POST activates nerve cells to inhibit pain.



UniTest® PC fits into a 40-well tray.



4X MAGNIFICATION

All test sites should show this impression on the skin with proper pressure.



NEG CNTL

UniTest® PC provides easy-to-read and highly reproducible histamine results.

FINALLY, CONTROLLED PAIN WITHOUT SMALL, HARD-TO-READ REACTIONS.

Made in USA

Decatur, Illinois • (800) 537-1336 • www.lincolndiagnostics.com
ISO 13485 certified • OSHA compliant • US and Foreign Patents Pending

L I N C O L N

Diagnosics, INC.